

Anti-Cystatin 9L2 (F87) CST9LP1 Antibody

Catalog Number: A19234-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Cystatin 9L2 (F87) CST9LP1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cystatin 9L2 (F87) CST9LP1 Antibody catalog # A19234-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q5W188

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 57-117 of Human Cystatin 9L2.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Cystatin 9L2 (F87) CST9LP1 Antibody (A19234-1) Images



Western blot (WB) analysis of Cystatin 9L2 (F87) polyclonal antibody at 1:500 dilution Lane1:HepG2 whole cell lysate(40ug) Lane2:HEK293T whole cell lysate(40ug) Lane3:C6 whole cell lysate(40ug) Lane4:3T3-L1 whole cell lysate(40ug)

[Anti-Cystatin 9L2 \(F87\) CST9LP1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.