

Anti-DHRS11 Antibody

Catalog Number: A30635

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-DHRS11 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-DHRS11 Antibody (Catalog# A30635). Tested in IHC, ICC, IF, WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6UWP2

Technical Details

Immunogen	A peptide derived from human DHRS11. Immunogen sequence location: 111-160
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Suggested Dilutions	WB 1:500-2000

	IHC 1:100 - 1:300. ELISA: 1:20000. ICC/IF 1:50-200
--	--

Anti-DHRS11 Antibody (A30635) Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using DHRS11 Antibody. The picture on the right is blocked with the synthesized peptide.

Anti-DHRS11 Antibody

For Research Use Only. Not for use in diagnostic procedures.