

Anti-DNA primase small subunit PRIM1 Antibody

Catalog Number: A30648

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-DNA primase small subunit PRIM1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DNA primase small subunit PRIM1 Antibody (Catalog# A30648). Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49642

Technical Details


Immunogen	A peptide derived from human PRIM1. Immunogen sequence location: 371-420
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Suggested Dilutions	WB 1:500 - 1:2000.

	IHC 1:100 - 1:300. ELISA: 1:10000. ICC/IF 1:50-200
--	--


Anti-DNA primase small subunit PRIM1 Antibody (A30648) Images




Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



Western Blot analysis of various cells using PRIM1 Polyclonal Antibody.



Western blot analysis of the lysates from K562 cells using PRIM1 antibody.



Western blot analysis of lysates from COLO and HeLa cells, using PRIM1 Antibody. The lane on the right is blocked with the synthesized peptide.

Anti-DNA primase small subunit PRIM1 Antibody

For Research Use Only. Not for use in diagnostic procedures.