

Anti-MARCH2 MARCH2; Antibody

Catalog Number: A30660

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-MARCH2 MARCH2; Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MARCH2 MARCH2; Antibody (Catalog# A30660). Tested in IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, ICC
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9P0N8


Technical Details

Immunogen	A peptide derived from human MARCH2. Immunogen sequence location: 191-240
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.


Suggested Dilutions

IHC 1:100 - 1:300.
ICC/IF 1:200 - 1:1000.
ELISA: 1:40000.
Not yet tested in other applications.

Anti-MARCH2 MARCH2; Antibody (A30660) Images



Immunofluorescence analysis of A549 cells, using MARCH2 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using MARCH2 Antibody. The picture on the right is blocked with the synthesized peptide.

[Anti-MARCH2 MARCH2; Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.