

Anti-MARCH3 MARCH3; Antibody

Catalog Number: A30661


Overview

Product Name	Anti-MARCH3 MARCH3; Antibody
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-MARCH3 MARCH3; Antibody (Catalog# A30661). Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat, Monkey.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q86UD3


Technical Details

Immunogen	A peptide derived from human MARCH3. Immunogen sequence location: 111-160
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Suggested Dilutions	WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:40000. ICC/IF 1:50-200


Anti-MARCH3 MARCH3; Antibody (A30661) Images




Western Blot analysis of K562 cells using MARCH3 Polyclonal Antibody



Western blot analysis of the lysates from RAW264.7 cells using MARCH3 antibody.



Immunohistochemistry analysis of paraffin-embedded human heart tissue, using MARCH3 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from K562 and COS7 cells, using MARCH3 Antibody. The lane on the right is blocked with the synthesized peptide.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MARCH3 MARCH3; Antibody

For Research Use Only. Not for use in diagnostic procedures.