

Anti-EIF3F/Eif3 Epsilon Antibody

Catalog Number: A30682

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-EIF3F/Eif3 Epsilon Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-EIF3F/Eif3 Epsilon Antibody (Catalog# A30682). Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O00303


Technical Details

Immunogen	A peptide derived from human EIF3F. Immunogen sequence location: 81-130
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.


Suggested Dilutions

WB 1:500 - 1:2000.
IHC 1:100 - 1:300.
ELISA: 1:40000.
ICC/IF 1:50-200


Anti-EIF3F/Eif3 Epsilon Antibody (A30682) Images



Western Blot analysis of various cells using eIF3epsilon Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using EIF3F Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HeLa cells, using EIF3F Antibody. The lane on the right is blocked with the synthesized peptide.

[Anti-EIF3F/Eif3 Epsilon Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.