

Anti-OR2AP1 Antibody

Catalog Number: A30868

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-OR2AP1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-OR2AP1 Antibody (Catalog# A30868). Tested in WB, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, ICC, WB
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NGE2

Technical Details


Immunogen	A peptide derived from human OR2AP1. Immunogen sequence location: 231 - 280
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Suggested Dilutions	WB 1:500 - 1:2000.

	ICC/IF 1:200 - 1:1000. ELISA: 1:5000. Not yet tested in other applications.
--	---


Anti-OR2AP1 Antibody (A30868) Images



Western Blot analysis of Jurkat cells using Olfactory receptor 2AP1 Polyclonal Antibody



Western blot analysis of the lysates from HT-29 cells using OR2AP1 antibody.



Western blot analysis of lysates from Jurkat cells, using OR2AP1 Antibody. The lane on the right is blocked with the synthesized peptide.

Anti-OR2AP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.