

Anti-Olfactory receptor 5AP2 OR5AP2 Antibody

Catalog Number: A30894

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Olfactory receptor 5AP2 OR5AP2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Olfactory receptor 5AP2 OR5AP2 Antibody (Catalog# A30894). Tested in WB, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, ICC, WB
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NGF4


Technical Details

Immunogen	A peptide derived from human OR5AP2. Immunogen sequence location: 267 - 316
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.


Suggested Dilutions

WB 1:500 - 1:2000.
ICC/IF 1:200 - 1:1000.
ELISA: 1:5000.
Not yet tested in other applications.

Anti-Olfactory receptor 5AP2 OR5AP2 Antibody (A30894) Images



Western Blot analysis of 3T3 cells using Olfactory receptor 5AP2 Polyclonal Antibody



Immunofluorescence analysis of HeLa cells, using OR5AP2 Antibody. The picture on the right is blocked with the synthesized peptide.

[Anti-Olfactory receptor 5AP2 OR5AP2 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.