

Anti-OR4A4/4A47 Antibody

Catalog Number: A30908


Overview

| | |
|----------------------|---|
| Product Name | Anti-OR4A4/4A47 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-OR4A4/4A47 Antibody (Catalog# A30908). Tested in WB, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat. |
| Application | ELISA, IF, ICC, WB |
| Clonality | Polyclonal |
| Formulation | PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q8NGN8 |


Technical Details

| | |
|---------------------|--|
| Immunogen | A peptide derived from human OR4A4. Immunogen sequence location: 250 - 299 |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure. |
| Purification | The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen. |
| Suggested Dilutions | WB 1:500 - 1:2000. ICC/IF 1:200 - 1:1000. ELISA: 1:20000. Not yet tested in other applications. |


Anti-OR4A4/4A47 Antibody (A30908) Images



Western blot analysis of the lysates from HT-29 cells using OR4A4/47 antibody.



Immunofluorescence analysis of A549 cells, using OR4A4/4A47 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from LOVO cells and HT29 cells, using OR4A4/4A47 Antibody. The lane on the right is blocked with the synthesized peptide.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-OR4A4/4A47 Antibody

For Research Use Only. Not for use in diagnostic procedures.