

Anti-ProSAPiP1/LZTS3 Antibody

Catalog Number: A31491

About LZTS3

Predicted to enable PDZ domain binding activity. Predicted to be involved in protein homooligomerization; regulation of dendritic spine morphogenesis; and regulation of postsynapse assembly. Predicted to be located in postsynaptic density. Predicted to be active in dendritic spine and glutamatergic synapse.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-ProSAPiP1/LZTS3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ProSAPiP1/LZTS3 Antibody catalog # A31491. Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O60299

Technical Details

Immunogen	E.coli-derived human ProSAPiP1/LZTS3 recombinant protein (Position: E330-T625).
Form	Liquid

Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ELISA, 1:100-1000

[Anti-ProSAPiP1/LZTS3 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.