

Anti-MB21D1 CGAS Antibody

Catalog Number: A31676


Overview

Product Name	Anti-MB21D1 CGAS Antibody
Reactive Species	Human
Description	Boster Bio Anti-MB21D1 CGAS Antibody catalog # A31676. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8N884

Technical Details

Immunogen	A synthetic peptide of human MB21D1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB :1:500~1:2000 IHC :1:20~1:100 ELISA: 1:5000~1:10000

Anti-MB21D1 CGAS Antibody (A31676) Images



Western blot analysis of extracts of various cell lines, using MB21D1 antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MB21D1 CGAS Antibody

For Research Use Only. Not for use in diagnostic procedures.