

Anti-KIAA0368/ECPAS Antibody

Catalog Number: A31777

About ECPAS

Ecm29 binds specifically to 26 S proteasomes. Ecm29 has been reported to stabilize proteasomes, to remodel proteasomes upon stress, to be involved in membrane-associated localization of proteasomes, and, more recently, to be involved in quality control or assembly of the proteasome. Ecm29 acts as a proteasome inhibitor.

Overview

Product Name	Anti-KIAA0368/ECPAS Antibody
Reactive Species	Human
Description	Boster Bio Anti-KIAA0368/ECPAS Antibody catalog # A31777. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q5VYK3

Technical Details

Immunogen	E.coli-derived human KIAA0368/ECPAS recombinant protein (Position: C1503-N1843).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KIAA0368/ECPAS Antibody

For Research Use Only. Not for use in diagnostic procedures.