

## Anti-PPP2R4 PTPA Antibody

Catalog Number: A31784

### Overview


Product Name	Anti-PPP2R4 PTPA Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-PPP2R4 PTPA Antibody catalog # A31784. Tested in ICC/IF applications. This antibody reacts with Human, Mouse.
Application	IF, ICC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15257

### Technical Details

Immunogen	Recombinant fusion protein of human Ser/thr-PP2A activator (PTPA)(NP_821067.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	ICC/IF: 1:50-1:200

## Anti-PPP2R4 PTPA Antibody (A31784) Images

---



Immunofluorescence analysis of U2OS cells using Ser/thr-PP2A activator antibody . Blue: DAPI for nuclear staining.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PPP2R4 PTPA Antibody

For Research Use Only. Not for use in diagnostic procedures.