

Anti-Brachyury TBXT Antibody

Catalog Number: A31984


Overview

Product Name	Anti-Brachyury TBXT Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-Brachyury TBXT Antibody catalog # A31984. Tested in WB,ICC/IF applications. This antibody reacts with Human,Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O15178


Technical Details

Immunogen	Recombinant Brachyury fusion protein containing a sequence corresponding to amino acids 226-435 of human Brachyury (NP_003172.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified.
Suggested Dilutions	WB: 1:500-1:1000 ICC/IF: 1:50-1:200

Anti-Brachyury TBXT Antibody (A31984) Images



Western blot analysis of H460 cells, using Brachyury Rabbit pAb antibody 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG antibody 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in Brachyury-BSA. Detection: ECL Basic Kit. Exposure time: 30s.



Immunofluorescence analysis of NCI-H460 cells using Brachyury Rabbit pAb antibody dilution of 1:200. Blue: DAPI for nuclear staining.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Brachyury TBXT Antibody

For Research Use Only. Not for use in diagnostic procedures.