

Anti-ADCK4 (G55) COQ8B Antibody

Catalog Number: A32263


Overview

Product Name	Anti-ADCK4 (G55) COQ8B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ADCK4 (G55) COQ8B Antibody catalog # A32263. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96D53

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 25-85 of Human ADCK4.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:5000-1:10000

Anti-ADCK4 (G55) COQ8B Antibody (A32263) Images



Western blot (WB) analysis of ADCK4 (G55) polyclonal antibody at 1:5000 dilution Lane1:SP2/0 whole cell lysate(30ug) Lane2:PC12 whole cell lysate(20ug) Lane3:Hela whole cell lysate(30ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ADCK4 (G55) COQ8B Antibody

For Research Use Only. Not for use in diagnostic procedures.