

Anti-ADCK3 (N323) COQ8A Antibody

Catalog Number: A32273

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-ADCK3 (N323) COQ8A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ADCK3 (N323) COQ8A Antibody catalog # A32273. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NI60

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 293-353 of Human ADCK3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 This affinity purified antibody has been tested for use in ELISA and western blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band 50-55 kDa in size corresponding to KSHV ORF57 protein by western blotting in the appropriate cell lysate or extract.

Anti-ADCK3 (N323) COQ8A Antibody (A32273) Images



Western blot (WB) analysis of ADCK3 (N323) polyclonal antibody at 1:500 dilution Lane1:CT26 whole cell lysate(40ug) Lane2:H9C2 whole cell lysate(40ug) Lane3:EC9706 whole cell lysate(40ug) Lane4:A2780 whole cell lysate(40ug) Lane5:Panc1 whole cell lysate(40ug)

[Anti-ADCK3 \(N323\) COQ8A Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.