

Anti-VISTA/VSIR Antibody

Catalog Number: A32320-1

About VSIR

VISTA, also known as platelet receptor Gi24, B7-H5, C10orf54, Dies1, PD-1H and SISP1, is a member of the immunoglobulin superfamily. It serves as not only as a positive regulator but also as a substrate of MT1-MMP, and contributes at least in part to tumor invasion.

Overview

Product Name	Anti-VISTA/VSIR Antibody
Reactive Species	Human
Description	Boster Bio Anti-VISTA/VSIR Antibody catalog # A32320-1. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9H7M9

Technical Details

Immunogen	E.coli-derived human VISTA/VSIR recombinant protein (Position: A36-A194).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunofluorescence, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-VISTA/VSIR Antibody

For Research Use Only. Not for use in diagnostic procedures.