



Isotype	IgG
Form	Liquid
Concentration	1 mg/mL
Purification	CCDC98 Antibody is affinity chromatography purified via peptide column.
Suggested Dilutions	CCDC98 antibody can be used for detection of CCDC98 by Western blot at 1 ug/mL. Antibody can also be used for immunohistochemistry starting at 5 ug/mL. Antibody validated: Western Blot in human samples and Immunohistochemistry in human samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.

Anti-CCDC98 ABRAXAS1 Antibody (A32338) Images



Western blot analysis of CCDC98 in human breast tissue lysate in (A) the absence and (B) presence of blocking peptide with CCDC98 antibody at 1 ug/mL.



Immunohistochemistry of CCDC98 in human breast tissue with CCDC98 antibody at 5 ug/mL.

[Anti-CCDC98 ABRAXAS1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.