

Anti-GRK2 (Phospho-Tyr86) Antibody

Catalog Number: A32388Y86


Overview

Product Name	Anti-GRK2 (Phospho-Tyr86) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GRK2 (Phospho-Tyr86) Antibody catalog # A32388Y86. Tested in IHC, IF, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P25098

Technical Details

Immunogen	Synthesized peptide derived from human GRK2 (Phospho-Tyr86)
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	IHC 1:50-200 WB 1:500-2000 IF 1:50-200

Anti-GRK2 (Phospho-Tyr86) Antibody (A32388Y86) Images



Immunohistochemical analysis of paraffin-embedded human small intestinal carcinoma tissue. 1, primary Antibody was diluted at 1:200 (4° overnight). 2, Sodium citrate pH 6.0 was used for antigen retrieval (>98°C, 20min). 3, Secondary antibody was diluted at 1:200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GRK2 (Phospho-Tyr86) Antibody

For Research Use Only. Not for use in diagnostic procedures.