

Anti-FYB FYB1 Antibody

Catalog Number: A32410-1

Overview

Product Name	Anti-FYB FYB1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-FYB FYB1 Antibody catalog # A32410-1. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O15117

Technical Details

Immunogen	Recombinant protein within C-terminal human FYB.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified.
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-FYB FYB1 Antibody (A32410-1) Images




Figure 2. Immunohistochemistry validation of FYB1 using Anti-FYB FYB1 Antibody (A32410-1).

Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-FYB antibody. The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 8.0-8.4) for 20 minutes. The tissues were blocked in 5% BSA for 30 minutes at room temperature




Figure 1. Western blotting validation for Anti-FYB FYB1 Antibody A32410-1

Western blot analysis of FYB on U937 cell lysate. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody was used at a 1:500 dilution in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG -HRP Secondary Antibody (HA1001) at 1:5

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FYB FYB1 Antibody