

## Anti-SKIV2L2 MTREX Antibody

Catalog Number: A32485


### Overview

Product Name	Anti-SKIV2L2 MTREX Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-SKIV2L2 MTREX Antibody catalog # A32485. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P42285


### Technical Details

Immunogen	Recombinant fusion protein of human SKIV2L2(NP_056175.3).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200 IP: 1:50-1:200


## Anti-SKIV2L2 MTREX Antibody (A32485) Images




Western blot analysis of extracts of various cell lines, using SKIV2L2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Immunohistochemistry of paraffin-embedded human breast cancer using SKIV2L2 antibody at dilution of 1:100. Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry of paraffin-embedded mouse testis using SKIV2L2 antibody at dilution of 1:100. Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunoprecipitation analysis of 300ug extracts of 293T cells using 3ug SKIV2L2 antibody. Western blot was performed from the immunoprecipitate using SKIV2L2 antibody at a dilution of 1:500.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SKIV2L2 MTREX Antibody

For Research Use Only. Not for use in diagnostic procedures.