

Anti-ND3 Mtnd3 Antibody

Catalog Number: A32721


Overview

Product Name	Anti-ND3 Mtnd3 Antibody
Reactive Species	Mouse
Description	Boster Bio Anti-ND3 Mtnd3 Antibody catalog # A32721. Tested in WB applications. This antibody reacts with Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P03899

Technical Details

Immunogen	A synthetic peptide corresponding to amino acids 21-34 of Mouse ND3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-ND3 Mtd3 Antibody (A32721) Images



Western blot (WB) analysis of ND3 polyclonal antibody at 1:1000 dilution Lane1:The muscle tissue lysate of Mouse (47ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ND3 Mtd3 Antibody

For Research Use Only. Not for use in diagnostic procedures.