

## Anti-APOBEC3D/F Antibody

Catalog Number: A32790


### Overview

Product Name	Anti-APOBEC3D/F Antibody
Reactive Species	Human, Monkey
Description	Boster Bio Anti-APOBEC3D/F Antibody catalog # A32790. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Monkey.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P60705


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human APOBEC3D/F. AA range:232-281
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:40000 IF 1:50-200

## Anti-APOBEC3D/F Antibody (A32790) Images



Immunohistochemistry analysis of paraffin-embedded human placenta tissue, using APOBEC3D/F Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of various cells using APOBEC3D/F Polyclonal Antibody

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-APOBEC3D/F Antibody

For Research Use Only. Not for use in diagnostic procedures.