

Anti-NADH2/Mtnd2 Antibody Picoband® HRP Conjugated

Catalog Number: A32839-HRP

About Mtnd2

Mitochondrially encoded NADH dehydrogenase 2 is protein that in humans is encoded by the mitochondrial gene MT-ND2 gene. The ND2 protein is a subunit of NADH dehydrogenase (ubiquinone), which is located in the mitochondrial inner membrane and is the largest of the five complexes of the electron transport chain. Variants of MT-ND2 are associated with mitochondrial encephalomyopathy, lactic acidosis, and stroke-like episodes (MELAS), Leigh's syndrome (LS), Leber's hereditary optic neuropathy (LHON) and increases in adult BMI.

Overview

Product Name	Anti-NADH2/Mtnd2 Antibody Picoband® HRP Conjugated
Reactive Species	Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P03893

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of mouse NADH2/Mtnd2, which shares 80% amino acid (aa) sequence identity with rat NADH2/Mtnd2.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NADH2/Mtnd2 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.