

Anti-ND1 Mtnd1 Antibody

Catalog Number: A32847


Overview

Product Name	Anti-ND1 Mtnd1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ND1 Mtnd1 Antibody catalog # A32847. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P03888

Technical Details

Immunogen	A synthetic peptide corresponding to amino acids 189-202 of Mouse ND1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-ND1 Mtdn1 Antibody (A32847) Images



Western blot (WB) analysis of ND1 polyclonal antibody at 1:1000 dilution Lane1:HEK293T whole cell lysate (57ug) Lane2:HepG-2 whole cell lysate(50ug) Lane3:H9C2 whole cell lysate(37ug) Lane4:NIH/3T3 whole cell lysate(39ug) Lane5:The muscle tissue lysate of Mouse (33ug) Lane6:The muscle tissue lysate of Rat(37ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ND1 Mtdn1 Antibody

For Research Use Only. Not for use in diagnostic procedures.