

Anti-LY6C Antibody Picoband® Biotin Conjugated

Catalog Number: A33313-Biotin

About Ly6c1/2

Ly6C is a murine cell-surface antigen expressed by plasma cells, subsets of myeloid cells and many T cells, including memory T cells.

Overview

Product Name	Anti-LY6C Antibody Picoband® Biotin Conjugated
Reactive Species	Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	POCW02

Technical Details

Immunogen	E.coli-derived mouse LY6C recombinant protein (Position: L27-G109).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LY6C Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.