

Anti-PLAC2 TINCR Antibody

Catalog Number: A33935

About TINCR

PLAC2 (placenta-specific protein 2) is a 151 amino acid protein that is encoded by a gene located on human chromosome 19p13. Based on its amino acid sequence, PLAC2 is predicted to be a multi-pass membrane protein. There are no reported studies or functions reported for this protein.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-PLAC2 TINCR Antibody
Reactive Species	Human
Description	Boster Bio Anti-PLAC2 TINCR Antibody (Catalog # A33935). Tested in ELISA, WB, ICC, IF applications. This antibody reacts with Human.
Application	ELISA, IF, ICC, WB
Clonality	Polyclonal
Formulation	PLAC2 Antibody is supplied in PBS containing 0.02% sodium azide.
Storage Instructions	PLAC2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. Avoid repeated freeze-thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Host	Rabbit
Uniprot ID	Q9NPR7

Technical Details


Immunogen	PLAC2 antibody was raised against a 17 amino acid synthetic peptide near the center of human PLAC2. The immunogen is located within amino acids 60 - 110 of PLAC2.
Predicted Reactive Species	Human
Isotype	IgG
Form	Liquid

Concentration	1 mg/mL
Purification	PLAC2 Antibody is affinity chromatography purified via peptide column.
Suggested Dilutions	PLAC2 antibody can be used for detection of PLAC2 by Western blot at 1 - 2 ug/mL. Antibody can also be used for immunocytochemistry starting at 10 ug/mL. For immunofluorescence start at 20 ug/mL. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.


Anti-PLAC2 TINCR Antibody (A33935) Images



Western blot analysis of PLAC2 in HeLa cell lysate with PLAC2 antibody at 1 ug/mL.



Immunocytochemistry of PLAC2 in HeLa cells with PLAC2 antibody at 10 ug/mL.



Immunofluorescence of PLAC2 in HeLa cells with PLAC2 antibody at 20 ug/mL.

[Anti-PLAC2 TINCR Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.