

Anti-MAGD4 MAGED4MAGED4B Antibody

Catalog Number: A33944

Overview

Product Name	Anti-MAGD4 MAGED4MAGED4B Antibody
Reactive Species	Human
Description	Boster Bio Anti-MAGD4 MAGED4MAGED4B Antibody catalog # A33944. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96JG8

Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 300-380
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	MAGD4 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MAGD4 MAGED4MAGED4B Antibody

For Research Use Only. Not for use in diagnostic procedures.