

## Anti-Casein Kinase epsilon CSNK1ELOC400927-CSNK1E Antibody

Catalog Number: A33958


### Overview

Product Name	Anti-Casein Kinase epsilon CSNK1ELOC400927-CSNK1E Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Casein Kinase epsilon CSNK1ELOC400927-CSNK1E Antibody catalog # A33958. Tested in ELISA, IF, IHC, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49674


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human CKI-epsilon. AA range:276-325
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:20000


## Anti-Casein Kinase epsilon CSNK1ELOC400927-CSNK1E Antibody (A33958) Images




Western Blot analysis of LOVO cells using Casein Kinase epsilon Polyclonal Antibody



Western blot analysis of lysates from LOVO cells, using CKI-epsilon Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human placenta tissue, using CKI-epsilon Antibody. The picture on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of HeLa cells, using CKI-epsilon Antibody. The picture on the right is blocked with the synthesized peptide.

**Submit a product review to [Biocompare.com](https://www.biocompare.com)**

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Casein Kinase epsilon CSNK1ELOC400927-CSNK1E Antibody

For Research Use Only. Not for use in diagnostic procedures.