

Anti-SARS-CoV-2 N/ORF9/Nucleocapsid Antibody

Catalog Number: A33996

About N

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is an enveloped, positive-sense, single-stranded RNA virus that causes coronavirus disease 2019 (COVID-19). Virus particles include the RNA genetic material and structural proteins needed for invasion of host cells. Once inside the cell the infecting RNA is used to encode structural proteins that make up virus particles, nonstructural proteins that direct virus assembly, transcription, replication and host control and accessory proteins whose function has not been determined.~ The structural proteins of SARS-CoV-2 include the envelope protein (E), spike or surface glycoprotein (S), membrane protein (M) and the nucleocapsid protein (N). The nucleocapsid phosphoprotein is a structural protein that binds to, protects the viral RNA genome and is involved in packaging the RNA into virus particles. The N protein has been suggested as an antiviral drug target.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-SARS-CoV-2 N/ORF9/Nucleocapsid Antibody
Reactive Species	Human
Description	Boster Bio Anti-SARS-CoV-2 N/ORF9/Nucleocapsid Antibody catalog # A33996. Tested in ELISA applications. This antibody reacts with Human.
Application	ELISA
Clonality	Polyclonal
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl and 0.2mg Na2HPO4.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P0DTC9

Technical Details

Immunogen	MSDNGPQNQRNAPRITFGGPSDESTGSNQNGERSGARSKQRRPQGLPNNTASWFTALTQHGKEDLK FPRGQGVPIINTNSSPDDQIGYYRRATRRIRGGDGKMKDLSPRWYFYLLGTGPEAGLPYGANKDGIW VATEGALNTPKDHIGTRNPANNAIIVLQLPQGTTLPKGFYAEGSRGGSQASSRSSRSRNSRNSTPG
-----------	---

	SSRGTSPARMAGNGGDAALALLLDRLNQLESKMSGKGQQQGGQTVTKKSAAEASKKPRQKRTAT KAYNVTQAFGRRGPEQTQGNFGDQELIRQGTDYKHWPQIAQFAPSASAFFGMSRIGMEVTPSGTW LTYTGAIKLDDKDPNFKDQVILLNKHIDAYKTFPPTEPKDKKKKKADETQALPQRQKKQQTVTLLPA ADLDDFSKQLQQSMSSADSTQA
Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	ELISA, 0.001-0.1ug/ml, Human

[Anti-SARS-CoV-2 N/ORF9/Nucleocapsid Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.