

## Anti-SARS-CoV-2 NSP4 Antibody

Catalog Number: A34000

### About rep

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is an enveloped, positive-sense, single-stranded RNA virus that causes coronavirus disease 2019 (COVID-19). Virus particles include the RNA genetic material and structural proteins needed for invasion of host cells. Once inside the cell the infecting RNA is used to encode structural proteins that make up virus particles, nonstructural proteins that direct virus assembly, transcription, replication and host control and accessory proteins whose function has not been determined.~ ORF1ab, the largest gene, contains overlapping open reading frames that encode polyproteins PP1ab and PP1a. The polyproteins are cleaved to yield 16 nonstructural proteins, NSP1-16. Production of the longer (PP1ab) or shorter protein (PP1a) depends on a -1 ribosomal frameshifting event. The proteins, based on similarity to other coronaviruses, include the papain-like proteinase protein (NSP3), 3C-like proteinase (NSP5), RNA-dependent RNA polymerase (NSP12, RdRp), helicase (NSP13, HEL), endoRNase (NSP15), 2'-O-Ribose-Methyltransferase (NSP16) and other nonstructural proteins. SARS-CoV-2 nonstructural proteins are responsible for viral transcription, replication, proteolytic processing, suppression of host immune responses and suppression of host gene expression. The RNA-dependent RNA polymerase is a target of antiviral therapies.

### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-SARS-CoV-2 NSP4 Antibody   |
| Reactive Species     | Human   |
| Description          | Boster Bio Anti-SARS-CoV-2 NSP4 Antibody catalog # A34000. Tested in ELISA applications. This antibody reacts with Human.   |
| Application          | ELISA   |
| Clonality            | Polyclonal  |
| Formulation          | Each vial contains 4mg Trehalose, 0.9mg NaCl and 0.2mg Na2HPO4.   |
| Storage Instructions | Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit  |
| Uniprot ID           | P0DTC1/P0DTD1   |

### Technical Details

|           |   |
|-----------|---|
| Immunogen | HVMSKHTDFSEIIGYKAIDGGVTRDIASDTDCFANKHADFDTWFSQRGGSYTNDKAggggDFLHFLP RVFSAVGNICYTPSKLIEYDFATSACVLAAECTIFKDasGKVPYCYDTNVLEGSVAYESLRPDTRggg gKHFYWFFSNYLKRRVVFNGVSFSTFEEAALCTFLLNKEMYLKLRSDVLLPLTQYNRYLALYNKYKYFS GAMDTTSYREAACCHLAKALNDFNSGSDVLYQPPQTSITSAVLQ |
| Isotype   | Rabbit IgG  |
| Form      | Lyophilized   |

|                     |   |
|---------------------|---|
| Concentration       | Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml. |
| Purification        | Immunogen affinity purified.  |
| Suggested Dilutions | ELISA, 0.001-0.1ug/ml, Human  |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SARS-CoV-2 NSP4 Antibody

For Research Use Only. Not for use in diagnostic procedures.