

## Anti-Zebrafish FOXO6A Antibody Picoband®

Catalog Number: AZA0A0R4IU78

### About FOXO6A

Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to act upstream of or within regulation of DNA-templated transcription. Predicted to be located in cytoplasm. Predicted to be active in nucleus. Is expressed in central nervous system; eye; hatching gland; and heart tube. Orthologous to human FOXO6 (forkhead box O6).

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



## Overview


|                      |                                                                                                                                                                                                                                                                                                                                                                                                                               |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Zebrafish FOXO6A Antibody Picoband®                                                                                                                                                                                                                                                                                                                                                                                      |
| Reactive Species     | Zebrafish                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Description          | Boster Bio Anti-Zebrafish FOXO6A Antibody Picoband® catalog # AZA0A0R4IU78. Tested in WB applications. This antibody reacts with Zebrafish. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance. |
| Application          | WB                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Formulation          | Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na <sub>2</sub> HPO <sub>4</sub> .                                                                                                                                                                                                                                                                                                                                     |
| Storage Instructions | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.                                                                                                                                                                                                                             |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Uniprot ID           | A0A0R4IU78                                                                                                                                                                                                                                                                                                                                                                                                                    |

## Technical Details

|           |                                                                          |
|-----------|--------------------------------------------------------------------------|
| Immunogen | E.coli-derived Zebrafish FOXO6A recombinant protein (Position: M1-E527). |
| Form      | Lyophilized                                                              |

|                     |                                                                           |
|---------------------|---------------------------------------------------------------------------|
| Concentration       | Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml. |
| Purification        | Immunogen affinity purified.                                              |
| Suggested Dilutions | Western blot, 0.25-0.5 ug/ml, Zebrafish                                   |

## Anti-Zebrafish FOXO6A Antibody Picoband® (AZA0A0R4IU78) Images



Western blot analysis of FOXO6A using anti-FOXO6A antibody (AZA0A0R4IU78). Electrophoresis was performed on a 10% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: whole male zebrafish tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-FOXO6A antigen affinity purified polyclonal antibody (AZA0A0R4IU78) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for FOXO6A at approximately 66 kDa. The expected band size for FOXO6A is at 66 kDa.

### Anti-Zebrafish FOXO6A Antibody

For Research Use Only. Not for use in diagnostic procedures.