

## Anti-Zebrafish TTC28 Antibody APC Conjugated

Catalog Number: AZA0A8M6Z8C1-APC

### About TTC28

Acts upstream of or within microtubule polymerization or depolymerization. Is expressed in blastomere and nervous system. Orthologous to human TTC28 (tetratricopeptide repeat domain 28).

### Overview

Product Name	Anti-Zebrafish TTC28 Antibody APC Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (IHC). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	A0A8M6Z8C1

### Technical Details

Immunogen	E.coli-derived zebrafish TTC28 recombinant protein (Position: Q1762-D1944).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish TTC28 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.