

Anti-Zebrafish IRF3 Antibody Picoband® Cy3 Conjugated

Catalog Number: AZA0A8M6Z9B8-Cy3

About IRF3

Enables transcription cis-regulatory region binding activity. Predicted to be involved in immune system process and regulation of transcription by RNA polymerase II. Predicted to act upstream of or within regulation of DNA-templated transcription. Predicted to be active in nucleus. Is expressed in head kidney. Human ortholog(s) of this gene implicated in brain disease. Orthologous to human IRF3 (interferon regulatory factor 3).

Overview

Product Name	Anti-Zebrafish IRF3 Antibody Picoband® Cy3 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	A0A8M6Z9B8

Technical Details

Immunogen	E.coli-derived Zebrafish IRF3 recombinant protein (Position: S92-L312).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish IRF3 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.