

APC Conjugated

Catalog Number: -APC

[description_before_attributes] !background!

About [gene_name]

[background]

!/background!

Overview

Product Name	APC Conjugated		
Reactive Species	[reactivity]		
!conjugate!	Conjugate	[conjugate]	
!/conjugate!	Application	Recommended applications are based on the parent unconjugated antibody (). Consult the data sheet for suitable applications according to their experimental needs.	
!clonality!	Clonality	[clonality]	
!/clonality!	!contents!	Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.05% BSA, and 0.05% sodium azide.
!/contents!	Storage Instructions		At -20°C for one year from date of manufacture.
!host!	Host	[host]	
!/host!	!uniprot_id!	Uniprot ID	[uniprot_id]
!/uniprot_id!			

Technical Details

!immunogen!	Immunogen	[immunogen]		
!/immunogen!	!predicted_reactivity!	Predicted Reactive Species	[predicted_reactivity]	
!/predicted_reactivity!	!cross_reactivity!	Cross Reactivity	[cross_reactivity]	
!/cross_reactivity!	!isotype!	Isotype	[isotype]	
!/isotype!	Formulation		[formulation]	
!purification!	Purification	[purification]		

!purification!	Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
	Suggested Dilutions	Optimal dilutions should be determined by end users.

[publication_count]

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



- APC

For Research Use Only. Not for use in diagnostic procedures.