

Anti-Zebrafish PKLR Antibody Picoband® Biotin Conjugated

Catalog Number: AZA0A8M9NZ14-Biotin

About PKLR

Predicted to enable pyruvate kinase activity. Acts upstream of or within liver development. Predicted to be active in cytoplasm. Is expressed in digestive system; female organism; heart; male organism; and muscle. Human ortholog(s) of this gene implicated in Gaucher's disease; congenital nonspherocytic hemolytic anemia; congenital nonspherocytic hemolytic anemia 2; malaria; and type 2 diabetes mellitus. Orthologous to human PKLR (pyruvate kinase L/R).

Overview

Product Name	Anti-Zebrafish PKLR Antibody Picoband® Biotin Conjugated
Reactive Species	Zebrafish
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	A0A8M9NZ14

Technical Details

Immunogen	E.coli-derived Zebrafish PKLR recombinant protein (Position: D31-D503).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.