

Anti-Zebrafish DUT Antibody Picoband® Biotin Conjugated

Catalog Number: AZA0A8M9P120-Biotin

About DUT

Predicted to enable dUTP diphosphatase activity and magnesium ion binding activity. Predicted to be involved in dUMP biosynthetic process and dUTP catabolic process. Predicted to act upstream of or within nucleotide metabolic process. Is expressed in digestive system; nervous system; neural plate; pectoral fin; and tail bud. Orthologous to human DUT (deoxyuridine triphosphatase).

Overview

Product Name	Anti-Zebrafish DUT Antibody Picoband® Biotin Conjugated
Reactive Species	Zebrafish
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	A0A8M9P120

Technical Details

Immunogen	E.coli-derived Zebrafish DUT recombinant protein (Position: K31-R157).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.