

## Anti-Zebrafish DUT Antibody Picoband® Fluoro647 Conjugated

Catalog Number: AZA0A8M9P120-Fluoro647

### About DUT

Predicted to enable dUTP diphosphatase activity and magnesium ion binding activity. Predicted to be involved in dUMP biosynthetic process and dUTP catabolic process. Predicted to act upstream of or within nucleotide metabolic process. Is expressed in digestive system; nervous system; neural plate; pectoral fin; and tail bud. Orthologous to human DUT (deoxyuridine triphosphatase).

### Overview

Product Name	Anti-Zebrafish DUT Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	A0A8M9P120

### Technical Details

Immunogen	E.coli-derived Zebrafish DUT recombinant protein (Position: K31-R157).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish DUT Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.