

Anti-Zebrafish TARDBPa/b Antibody Picoband® Fluoro488 Conjugated

Catalog Number: AZA0A8M9PSZ1-Fluoro488

About TARDBPa/b

Predicted to enable RNA binding activity. Acts upstream of or within several processes, including axon extension; neuromuscular junction development; and skeletal muscle tissue development. Predicted to be located in mitochondrion and nucleus. Predicted to be active in chromatin and nucleoplasm. Is expressed in brain and fin. Human ortholog(s) of this gene implicated in Parkinson's disease; amyotrophic lateral sclerosis; amyotrophic lateral sclerosis type 10; and motor neuron disease. Orthologous to human TARDBP (TAR DNA binding protein).

Overview

Product Name	Anti-Zebrafish TARDBPa/b Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	A0A8M9PSZ1/Q802C7

Technical Details

Immunogen	E.coli-derived zebrafish TARDBPa/b recombinant protein (Position: T2-H265).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish TARDBPa/b Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.