

Anti-Zebrafish MMP19 Antibody Fluoro594 Conjugated

Catalog Number: AZA0AB32TZ88-Fluoro594

About MMP19

Orthologous to human MMP19 (matrix metalloproteinase 19).

Overview

Product Name	Anti-Zebrafish MMP19 Antibody Fluoro594 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (IHC). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	A0AB32TZ88

Technical Details

Immunogen	E.coli-derived Zebrafish MMP19 recombinant protein (Position: A237-K396).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.