

Anti-Zebrafish FOXD1 Antibody HRP Conjugated

Catalog Number: AZA2CE80-HRP

About FOXD1

Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in anatomical structure morphogenesis; cell differentiation; and regulation of transcription by RNA polymerase II. Predicted to act upstream of or within regulation of DNA-templated transcription. Predicted to be located in nucleus. Is expressed in immature eye; nervous system; neural crest; pharyngeal arch; and presumptive diencephalon. Orthologous to human FOXD1 (forkhead box D1).

Overview

Product Name	Anti-Zebrafish FOXD1 Antibody HRP Conjugated
Reactive Species	Zebrafish
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	A2CE80

Technical Details

Immunogen	E.coli-derived Zebrafish FOXD1 recombinant protein (Position: D29-C343).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish FOXD1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.