

Anti-Zebrafish IHHA Antibody Picoband® HRP Conjugated

Catalog Number: AZA3KNS3-HRP

About IHHA

Predicted to enable calcium ion binding activity and patched binding activity. Acts upstream of or within several processes, including enteric nervous system development; intramembranous ossification; and swim bladder morphogenesis. Predicted to be located in several cellular components, including Golgi membrane; endoplasmic reticulum membrane; and extracellular region. Predicted to be active in extracellular space. Is expressed in several structures, including caudal fin skeleton; cranium; digestive system; pneumatic duct; and swim bladder bud. Human ortholog(s) of this gene implicated in acrocapitofemoral dysplasia; brachydactyly type A1; retinopathy of prematurity; and syndactyly type 1. Orthologous to human IHH (Indian hedgehog signaling molecule).

Overview

Product Name	Anti-Zebrafish IHHA Antibody Picoband® HRP Conjugated
Reactive Species	Zebrafish
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	A3KNS3

Technical Details

Immunogen	E.coli-derived Zebrafish IHHA recombinant protein (Position: A203-R413).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish IHHA Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.