

## Anti-Zebrafish GNRHR4 Antibody Fluoro647 Conjugated

Catalog Number: AZA3QJZ0-Fluoro647

### About GNRHR4

Enables gonadotropin-releasing hormone receptor activity. Acts upstream of or within G protein-coupled receptor signaling pathway. Predicted to be located in membrane. Predicted to be active in plasma membrane. Is expressed in several structures, including female organism; gonad; gut; male organism; and nervous system. Orthologous to human GNRHR2 (gonadotropin releasing hormone receptor 2 (pseudogene)).

### Overview

Product Name	Anti-Zebrafish GNRHR4 Antibody Fluoro647 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (IHC). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	A3QJZ0

### Technical Details

Immunogen	E.coli-derived Zebrafish GNRHR4 recombinant protein (Position: T7-E406).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish GNRHR4 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.