

Anti-Zebrafish PHYHIPL/b Antibody HRP Conjugated

Catalog Number: AZA4QNW7-HRP

About PHYHIPL/b

Phytanoyl-CoA hydroxylase-interacting protein-like is an enzyme that in humans is encoded by the PHYHIPL gene. PHYHIPL (phytanoyl-CoA 2-hydroxylase interacting protein-like), also known as phytanoyl-CoA hydroxylase-interacting protein-like, is a 376 amino acid protein that contains one fibronectin type-III domain and belongs to the PHYHIP family. Conserved in chimpanzee, canine, mouse, rat, chicken, zebrafish and *Caenorhabditis elegans*, PHYHIPL exists as three alternatively spliced isoforms. PHYHIPL is a down-regulated target of IRX1, a homeobox tumor suppressor gene linked to gastric carcinoma. PHYHIPL may also play a role in the development of the central system. The gene that encodes PHYHIPL maps to human chromosome 10q21.1.

Overview

Product Name	Anti-Zebrafish PHYHIPL/b Antibody HRP Conjugated
Reactive Species	Zebrafish
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	A4QNW7/A0A8M2B4A8

Technical Details

Immunogen	E.coli-derived zebrafish PHYHIPL/b recombinant protein (Position: M1-D337).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish PHYHIPL/b Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.