

Anti-Zebrafish DOCK1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: AZA8E1V6-Fluoro488

About DOCK1

Predicted to enable guanyl-nucleotide exchange factor activity and small GTPase binding activity. Acts upstream of or within myelination in peripheral nervous system and myoblast fusion. Predicted to be active in cytoplasm and plasma membrane. Orthologous to human DOCK1 (dedicator of cytokinesis 1).

Overview

Product Name	Anti-Zebrafish DOCK1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	A8E1V6

Technical Details

Immunogen	E.coli-derived Zebrafish DOCK1 recombinant protein (Position: S264-E1722).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish DOCK1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.