

Anti-Zebrafish VPS53 Antibody

Catalog Number: AZE7FAG0

About VPS53

Predicted to be involved in retrograde transport, endosome to Golgi. Predicted to be located in Golgi apparatus and endosome membrane. Predicted to be part of GARP complex. Human ortholog(s) of this gene implicated in pontocerebellar hypoplasia type 2E. Orthologous to human VPS53 (VPS53 subunit of GARP complex).


Overview

| | |
|----------------------|---|
| Product Name | Anti-Zebrafish VPS53 Antibody |
| Reactive Species | Zebrafish |
| Description | Boster Bio Anti-Zebrafish VPS53 Antibody catalog # AZE7FAG0. Tested in IHC applications. This antibody reacts with Zebrafish. |
| Application | IHC |
| Clonality | Polyclonal |
| Formulation | Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ . |
| Storage Instructions | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | E7FAG0 |


Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived zebrafish VPS53 recombinant protein (Position: E11-L695). |
| Form | Lyophilized |
| Concentration | Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml. |
| Purification | Immunogen affinity purified. |
| Suggested Dilutions | Immunohistochemistry(Paraffin-embedded Section), 2-5 ug/ml, Zebrafish Immunofluorescence, 2 ug/ml, Zebrafish |


Anti-Zebrafish VPS53 Antibody (AZE7FAG0) Images



IHC analysis of VPS53 using anti-VPS53 antibody (AZE7FAG0). VPS53 was detected in a paraffin-embedded section of zebrafish colon tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 2 ug/ml rabbit anti-VPS53 Antibody (AZE7FAG0) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.



IHC analysis of VPS53 using anti-VPS53 antibody (AZE7FAG0). VPS53 was detected in a paraffin-embedded section of zebrafish liver tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 2 ug/ml rabbit anti-VPS53 Antibody (AZE7FAG0) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.



IHC analysis of VPS53 using anti-VPS53 antibody (AZE7FAG0). VPS53 was detected in a paraffin-embedded section of zebrafish kidney tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 2 ug/ml rabbit anti-VPS53 Antibody (AZE7FAG0) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish VPS53 Antibody

For Research Use Only. Not for use in diagnostic procedures.