

Anti-Zebrafish SUB1a/b Antibody Picoband® Fluoro594 Conjugated

Catalog Number: AZQ504B1-Fluoro594

About SUB1a/b

Predicted to enable transcription coactivator activity. Predicted to act upstream of or within positive regulation of transcription initiation by RNA polymerase II. Predicted to be part of transcription regulator complex. Predicted to be active in nucleus. Is expressed in central nucleus of ventral telencephalon. Orthologous to human SUB1 (SUB1 regulator of transcription).

Overview

Product Name	Anti-Zebrafish SUB1a/b Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q504B1/Q6DGP6

Technical Details

Immunogen	E.coli-derived zebrafish SUB1a/b recombinant protein (Position: N58-T123).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish SUB1a/b Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.