

## Anti-Zebrafish CD326/EPCAM Antibody Picoband® Fluoro647 Conjugated

Catalog Number: AZQ568H0-Fluoro647

### Overview

Product Name	Anti-Zebrafish CD326/EPCAM Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q568H0

### Technical Details

Immunogen	E.coli-derived Zebrafish CD326/EPCAM recombinant protein (Position: S43-S303).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish CD326/EPCAM Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.