

## Anti-Zebrafish GIPC1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: AZQ6B337-Fluoro488

### About GIPC1

GIPC PDZ domain containing family, member 1 (GIPC1) is a protein that in humans is encoded by the GIPC1 gene. GIPC1 is a scaffolding protein that regulates cell surface receptor expression and trafficking.

### Overview

Product Name	Anti-Zebrafish GIPC1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Zebrafish
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6B337

### Technical Details

Immunogen	E.coli-derived zebrafish GIPC1 recombinant protein (Position: H70-V335).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.