

Anti-Zebrafish SUZ12a/b Antibody APC Conjugated

Catalog Number: AZQ6DC03-APC

About SUZ12a/b

Predicted to enable chromatin DNA binding activity. Predicted to act upstream of or within chromatin organization. Predicted to be part of ESC/E(Z) complex and RSC-type complex. Predicted to be active in nucleus. Is expressed in gonad; midbrain hindbrain boundary; pharyngeal arch; and visual system. Human ortholog(s) of this gene implicated in renal Wilms' tumor. Orthologous to human SUZ12 (SUZ12 polycomb repressive complex 2 subunit).

Overview

Product Name	Anti-Zebrafish SUZ12a/b Antibody APC Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (IHC). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6DC03/B0R1D5

Technical Details

Immunogen	E.coli-derived zebrafish SUZ12a/b recombinant protein (Position: K38-K448).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish SUZ12a/b Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.